

Cell name: MDCKII x ABCB1-Neo
Cat. No. M-0250

Product description


MDCKII x ABCB1-Neo was established by retroviral transduction of human ABCB1 (also known as P-glycoprotein or MDR1) cDNA into MCDKII cells and flow cytometry sorting of the calcein-extruding population.

Transgene construct information

Expression cassette type

Single-promoter, bicistronic

Expression cassette map



Expression cassette features

Element	Type	Species	RefSeq	Mutation / Discrepancy
SFFV	promoter, viral, constitutive	-	-	-
Splice donor site	Splice donor site	-	-	-
ABCB1	CDS	Homo sapiens	AF016535	none / none
Splice acceptor site	Splice acceptor site	-	-	-
Neomycin resistance gene	CDS	-	-	none / none

Transgene protein information

The transmembrane protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by this gene is an ATP-dependent drug efflux pump for xenobiotic compounds with broad substrate specificity. It is responsible for decreased drug accumulation in multidrug-resistant cells and often mediates the development of resistance to anticancer drugs. This protein also functions as a transporter in the blood-brain barrier. [Source: NCBI Gene, <http://www.ncbi.nlm.nih.gov/gene/5243>]

Creative Cell Ltd.

Homepage: www.creativecell.hu Email: info@creativecell.hu Phone: +36 (30) 41-44446

Cell culture characteristics

Source

Dog (*Canis familiaris*) kidney

Morphology

Epithelial

Growth properties

Adherent

Culture conditions

Culture cells at 37°C in humidified atmosphere with 5% CO₂. The base medium for this cell line is D-MEM with 4.5 g/l glucose. To make the complete growth medium, add fetal bovine serum to a final concentration of 10%.

Subculturing

Remove culture medium and detach cells by treating with 0.04 mL/cm² of 0.25% trypsin / 1 mM EDTA solution for 10 minutes at 37°C. Add 5x volume of complete medium to neutralize trypsin-EDTA, pellet cells at 250 x g, discard supernatant, and resuspend in culture medium. Subcultivate in a ratio of 1:4-1:6 two times a week. Recommended plating density: 3x10⁴ cells/cm².

Preservation

Freeze Medium: Complete growth medium with 10% DMSO.

Storage Temperature: Liquid nitrogen vapour phase.


Population doubling time

~24 hours

Sterility testing

Mycoplasma: negative

Validation results



A: Dotted line: MDCKII x ABCB1-Neo cells incubated with isotype control.

Red line: MDCKII x ABCB1-Neo cells incubated with anti-Human ABCB1 antibody.

Indirect immunofluorescent labeling with anti-Human ABCB1 antibody (clone MRK16, Alexis Biochemicals / Enzo Life Sciences) or isotype control + PE labeled anti-Mouse IgG secondary antibody.

Detection by flow cytometry.

B: Red line: MDCKII x ABCB1-Neo cells incubated with calcein and verapamil.

Dotted line: MDCKII x ABCB1-Neo cells incubated with calcein.

Calcein uptake was detected by flow cytometry at 485 nm in the presence or absence of verapamil (inhibitor of ABCB1).

Creative Cell Ltd.

Homepage: www.creativecell.hu Email: info@creativecell.hu Phone: +36 (30) 41-44446